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3. SCRIPT DE BASE DE DATOS (schema.sql)

-- 1) Tabla sessions:

```
CREATE TABLE IF NOT EXISTS sessions (  
  id SERIAL PRIMARY KEY,  
  filename TEXT NOT NULL,  
  loaded_at TIMESTAMP NOT NULL,  
  scan_rate DOUBLE PRECISION,  
  start_potential DOUBLE PRECISION,  
  end_potential DOUBLE PRECISION,  
  software_version TEXT  
);
```

-- 2) Tabla measurements:

```
CREATE TABLE IF NOT EXISTS measurements (
```

```

id SERIAL PRIMARY KEY,
session_id INTEGER NOT NULL REFERENCES sessions(id) ON DELETE
CASCADE,
title TEXT NOT NULL,
timestamp TIMESTAMP NOT NULL,
device_serial TEXT NOT NULL,
curve_count INTEGER NOT NULL,
pca_data DOUBLE PRECISION[] NOT NULL,
ppm_estimations DOUBLE PRECISION[] NOT NULL
);
-- 3) Tabla curves:
CREATE TABLE IF NOT EXISTS curves (
id SERIAL PRIMARY KEY,
measurement_id INTEGER NOT NULL REFERENCES measurements(id) ON
DELETE CASCADE,
curve_index INTEGER NOT NULL,
num_points INTEGER NOT NULL
);
-- 4) Tabla points:
CREATE TABLE IF NOT EXISTS points (
id SERIAL PRIMARY KEY,
curve_id INTEGER NOT NULL REFERENCES curves(id) ON DELETE CASCADE,
potential DOUBLE PRECISION NOT NULL,
current DOUBLE PRECISION NOT NULL
);
-- Índices:
CREATE INDEX IF NOT EXISTS idx_measurements_session ON
measurements(session_id);
CREATE INDEX IF NOT EXISTS idx_curves_measurement ON
curves(measurement_id);
CREATE INDEX IF NOT EXISTS idx_points_curve ON points(curve_id);

```

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Bloque 2.4 on_session_select(): lee de sessions scan_rate, start_potential, end_potential, software_version

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